INNOVATIVE TOOLS FOR

*cytogenetics*

*Amnio*® **Plus**  
*med*

EK AMG-200  
EK AMG- 600

INSTRUCTION  
FOR USE

CE

IVD

**Euro**  **clone**®

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## **Amniocyte Culture "In situ" Protocol**

This protocol provides a guide for human amniocyte culture using **Amniomed® Plus**. It can be used to replace either part or the whole of existing protocols for amniocyte culture at the user's discretion.

### **"In situ" Cell Culture**

- a) Transfer the amniotic fluid sample (about 20 ml ) in 2 sterile centrifuge tubes and spin at 1200 rpm for 10 minutes
- b) Decant the supernatant for Alpha-fetoprotein analysis, leaving about 1 ml in each tube.

*Please Note: Thaw the **Amniomed® Plus** in a  $37\pm 2$  °C waterbath and mix well by swirling prior to use. Addition of L-Glutamine, antibiotics and serum are not necessary since these components are already present.*

- c) Prepare 4- 6 sterile Petri dishes (35 mm diameter) containing round glass slide (AMNIODISH cod. EK AMN-240)
- d) Distribute cell suspension into Petri dishes.
- e) Add 2 ml of **Amniomed® Plus** to each Petri dish.
- f) Put the cultures in a  $37\pm 2$  °C at 5% CO<sub>2</sub> incubator.
- g) On the 5<sup>th</sup> day observe carefully the cell growing condition and then replace the medium with 2 ml of fresh one.
- h) One or two days later observe the cultures and if the cell colony growth is adequate add Colcemid solution at the final concentration of 0.05µg/ml.

Maintain in incubator for 2-4 additional hours before chromosome preparation.

### **Flask method culture**

(Use the same procedure as for the above "in situ" culture with the following adaptation)

- a) Distribute cell suspension into 2-4 flasks.
- b) Add to each flask 4-5 ml of **Amniomed® Plus**.

c) Place cap loosely on the flask and put the cultures in a  $37\pm 2$  °C at 5% CO<sub>2</sub> incubator.

### ***Chromosome preparation***

a) Decant the medium completely by tilting the plate and suctioning from the corner.

b) Add 3.0 ml of hypotonic solution in each Petri dish and maintain at room temperature for 10 minutes

c) Add 0.5 ml of fixative mixture (3:1 methanol – acetic acid) directly to the hypotonic solution and leave for 5 minutes.

d) Suction off the supernatant and immediately add 3.0 ml of fresh fixative.

e) Repeat twice to remove carefully any trace of residual water present in the Petri dishes. (The duration of these washings has no influence on final result).

f) Remove coverslip from the dish and immediately put it to dry in OPTICHROME device (Optichrome cod. EK AMH-950) at right condition of temperature and relative humidity.

### **Chorionic Villus Sampling (CVS) Culture Protocol**

a) Transfer the specimen from the transport tube to a 60 mm Petri dish containing 5 ml of RPMI 1640

b) Wash Villi with fresh medium to remove the blood cells

c) Using the inverted microscope, carefully dissect and remove any remaining clots or decidua from Chorionic Villi.

d) Wash carefully Villi in Hank's Balanced Salt Solution without Ca<sup>++</sup> and Mg<sup>++</sup>.

### **Long-term cultures**

a) Transfer the Villi in a sterile centrifuge tube containing 2 ml of Pronase E (Meerk, 4.000.000 PU/g) \* and incubate for 4-6 minutes at room temperature, gently shaking from time to time.

b) Add 3 - 5 ml of cold HBSS (without Ca<sup>++</sup> and Mg<sup>++</sup>) to stop enzyme action.

c) Centrifuge for 5 minutes at 1500 rpm and discard the supernatant

- d) Add 2 ml of sterile Collagenase type II (Sigma, 1 mg/ml) and incubate at 37 °C for 10 minutes
- e) Add 3 - 5 ml of cold HBSS to stop the enzyme action.
- f) Centrifuge for 5 minutes at 1500 rpm and discard the supernatant
- g) Add 2 ml of **Amniomed® Plus** and resuspend the pellet.
- h) Prepare 2 to 6 Petri dishes (Amniodish) depending on the volume of cell pellet. Add 2 ml of **Amniomed® Plus** medium to each Petri dish and distribute the cell suspension.
- i) Incubate at 37 °C in a 5% CO<sub>2</sub> incubator.
- l) After 4 days of culture, check for progress of growth using an inverted microscope
- m) When the areas of growth show a number of cells in mitosis, add one drop of Colcemid to each Petri dish (at a final concentration of 0.05 µg/ml ) and keep in incubator for additional 4 - 6 hours.

### **Chromosome preparation**

Follow the same procedure described in "Chromosome preparation" for Amniotic Fluid cell culture.

\* Dissolve 25 mg of Pronase E in 25 ml of HBSS without Ca<sup>++</sup> and Mg<sup>++</sup>. Filter with a 0.2 µm membrane. Aliquot and put in deep freezer (≤ - 18°C).

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Euroclone S.p.A. - Life Sciences Division