INNOVATIVE TOOLS FOR

cytogenetics

*Chromosome
kit and medium M*

EK AMT-200M
EK AMT-500M
EK AMTB-100M
EK AMTB-500M

**INSTRUCTION
FOR USE**

CE

IVD

Euro  **clone**[®]

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CHROMOSOME KIT M : INSTRUCTIONS FOR USE

1. Aseptically add up to 1 ml of bone marrow to the tube (the volume is determined considering the cell density of the sample); tighten the cap and mix
2. Incubate horizontally at $37 \pm 2^\circ \text{C}$ for 12, 24 or 48 hours (CO_2 not required)
3. Add 4-5 drops of Colcemid and incubate (time must be decided empirically)
4. Centrifuge at 2000 rpm for 3-4 minutes
5. Discard supernatant, leaving about 0.5 ml of pellet material
6. Resuspend vigorously (vortexing is recommended)
7. Add 5 ml of hypotonic solution KCl 0.075M* - Wait 10 minutes at RT (Room Temperature).
8. Repeat steps 4), 5), 6)
9. Add 5 ml of modified Ibraimov's solution*[^] (92 ml H_2O +5 ml acetic acid+3 ml methanol)
10. Repeat steps 4), 5), 6)
11. Add 5 ml of fresh Carnoy's fixative solution* - Wait 10 minutes at RT (3:1 methanol-glacial acetic acid)
12. Repeat steps 4), 5)
13. Add 5 ml of Carnoy's fixative solution*
14. Repeat step 4) and discard supernatant almost completely
15. Resuspend the pellet with a few drops of freshly prepared fixative
16. Drop onto clean, wet slides under climatically controlled conditions^o

Notes:

* Dispensing must be done energetically in order to ensure better uniformity of the cell suspension

[^] An homogeneous brownish colour and foam formation will indicate that this step has been performed correctly

^o The following procedure is strongly suggested to obtain optimal conditions:

- a) allow one drop to evaporate and then observe the slide using a phase contrast microscope; if the chromosomes appear well spread, with sharp, clear edges and no cytoplasmic residues, then proceed with the final preparation dispensing three drops on the slide;
- b) should you see metaphases with lack of chromosomes ("broken" metaphases) this is due to an excess of humidity: it is advisable to place your slides under a normal lamp turned on, in order to reduce humidity;
- c) if, viceversa, the chromosomes are dense and retaining cytoplasm all around, this is because the environment is too dry: you can obviate simply by means of an humidity source, such as a water bath.

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Euroclone S.p.A. - Life Sciences Division