
Chromoprobe Multiprobe[®] - System OctoChrome[™]

REF PMP803/ PMP804/PMP802

ENGLISH/FRANÇAIS/ITALIANO/
DEUTSCH/ESPAÑOL/DANSK

ENGLISH

Contents

Introduction	4
Material provided	5
Warning and Precautions	5
Storage and Handling	5
Ancillary Products	5
Materials Necessary but not Supplied	5
Optimal Microscope and Filter Set-up	6
Sample Preparation	7
Chromoprobe Multiprobe Protocol	7
Slide preparation	7
Preparation of OctoChrome device and template slide	9
Positioning of the template over the OctoChrome device	9
Instructions for use of the Cytocell slide surface thermometer	10
Denaturation	10
Hybridisation	10
Post-hybridisation stringent washes	11
Mounting and visualisation of results	11
Counterstaining with PI or DAPI	11
Removal of PI and re-staining with DAPI	12
Stability of Finished Slides	12
Procedural Recommendations	12
Expected Results	12
Probe Specifications	14
References	14
FRANÇAIS	17
ITALIANO	31
DEUTSCH	45
ESPAÑOL	61
DANSK	75


Introduction

Fluorescence *In Situ* Hybridisation (FISH) is a technique that allows DNA sequences to be detected on metaphase chromosomes or in interphase nuclei of fixed cultured or uncultured cytogenetic samples. The technique uses DNA probes that hybridise to entire chromosomes or single unique sequences, and serves as a powerful adjunct to classic cytogenetics. Target DNA, after fixation, is treated with heat and formamide to denature the double-stranded DNA, rendering it single-stranded. The target DNA is thus available for annealing to a similarly denatured, single-stranded, fluorescently labelled DNA probe which has a complementary sequence. Following hybridisation, unbound and non-specifically bound DNA probe is removed by a series of stringent washes and the DNA counterstained for visualisation. Fluorescence microscopy then allows the visualisation of the hybridised probe on the target material.

The OctoChrome combines the principles of the Chromoprobe Multiprobe[®] System with multiple colour fluorescence labelling to provide a FISH technology to aid cytogenetic analysis, making it quick, simple and cost effective.

The Multiprobe-OctoChrome device is divided into 8 square areas, each containing three different whole chromosome painting probes each directly labelled in a different colour : green (FITC spectrum), red (Texas Red spectrum), and blue (Coumarin filter). The probes are reversibly bound to a glass device, a method exclusive to Cytocell, thereby simplifying fluorescence *in situ* hybridisation by removing the need to prepare the probes. Denaturation of probe and target DNA occurs simultaneously under the device once heated.

The FISH protocol is further simplified with the simultaneous denaturation of both probe and target DNA and, after overnight hybridisation, the use of rapid formamide-free stringency washes. The directly labelled painting probes remove the need for lengthy amplification steps.

The Multiprobe-OctoChrome is intended for FISH to metaphase chromosomes from fixed cultured peripheral blood cells.

Material Provided

Each kit contains the following reagents, which are sufficient for either 2 (Cat. No. PMP802), 5 (Cat. No. PMP804) or 10 (Cat. No. PMP803) patient samples:

- 2, 5 or 10 Chromoprobe Multiprobe - OctoChrome devices coated with directly labelled painting probes.
Amount of painting probe labelled in red: 2.5 to 5 ng per square.
Amount of painting probe labelled in green: 7.5 to 12.5 ng per square.
Amount of painting probe labelled in blue: 15 to 25 ng per square.
- 4, 7 or 12 Glass slides printed with a special template
- 500µl Hybridisation Solution A: Formamide, Dextran Sulphate, SSC
- 500µl Counterstain Solution: PI (ES: 0.125µg/ml PI (Propidium Iodide)), Antifade
- 500µl Counterstain Solution: DAPI (ES: 0.125µg/ml DAPI (4,6-diamidino-2-phenylindole)), Antifade
- 1 Cytocell Slide Surface Thermometer
- 1 Cytocell Chromoprobe Multiprobe Hybridisation Chamber

Warnings and precautions

1. For *in vitro* diagnostic use. For professional use only.
2. Wear gloves when handling DNA probes and DAPI counterstain.
3. The hybridisation solution contains formamide, which is a toxic. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
4. DAPI and PI are potential carcinogens. Handle with care; wear gloves and a lab coat. Upon disposal, flush with a large volume of water.
5. All hazardous materials should be disposed of according to your institution's guidelines for hazardous waste disposal.

Storage and Handling

The Chromoprobe Multiprobe System kit should be stored at 2-8°C until the expiry date indicated on the kit label. **Do not freeze.** The counterstain vial must be stored in the dark.

Ancillary Products

- 5 x ST Buffer - 100 ml Cat. No. PCA 004
- 20 x SSC Buffer Cat. No. PCA 003
- Cytocell Hotplate Cat. No. PCN 001
- Fluorescent filter sets are available from Cytocell upon request.

Materials Necessary but not Supplied

Equipment

- a) Hotplate with accurate temperature control up to 80°C
- b) 37°C incubator
- c) Variable volume micropipettes range 1µl - 200µl
- d) 37°C water bath (without stirrer)

- e) Microcentrifuge tubes (0.5 ml)
- f) Waterbath with accurate temperature control at 72°C
- g) Fluorescence microscope (Please see Optimal microscope and filter set up)
- h) Plastic or glass coplin jars
- i) Centrifuge
- j) Forceps
- k) Fluorescence grade microscope lens immersion oil
- l) Fluorescence grade glass coverslips (24 x 50mm)

Solutions

- a) **Carnoy's Fixative**
3:1 methanol : glacial acetic acid
- b) **Ethanol Series**
70%, 85%, 100% (absolute) ethanol
- c) **20 x SSC Stock**
175.3g Sodium chloride, 88.2g Sodium citrate per litre adjust to pH 7.0 with conc. HCl. (4 x SSC: dilute stock 1:5)(2 x SSC: dilute stock 1:10). This solution may be stored for up to 1 year at 4°C
- d) **ST Buffer (1x)**
4 x SSC plus 0.05% (v/v) Tween 20
- e) **0.4 x SSC**
Mix 2ml 20 x SSC with 98mls distilled water.
- f) **2 x SSC, 0.05% Tween 20**
Mix 10ml 20 x SSC with 90mls distilled water, add 50µl Tween 20
- g) **100% Methanol**

Optimal Microscope and Filter Set up

To ensure optimal results with this product it is essential that the microscope used is set up as follows:

General Microscope Set up	
100 Watt Mercury bulb	x 100 ; x 60 Plan-Apochromat Objectives

Fluorophore or counterstain visible

Filter Specifications	Green	Red	Blue with DAPI	Blue with PI	DAPI	PI
DAPI/ FITC/ Texas Red Triple Filter	v	v	v	v	v	v
FITC (narrow bandpass) Single Bandpass Filter	v					
Texas Red Single Bandpass Filter		v				v
Aqua Single Bandpass Filter			v	v		

Sample preparation

The Chromoprobe Multiprobe System is designed for use on cultured peripheral blood cells fixed in Carnoy's fixative which should be prepared according to the laboratory or institution guidelines. (See also Solutions, above)

Prepare blood metaphase spreads on CytoCELL Chromoprobe Multiprobe template slides according to CytoCELL protocol below. Baking or otherwise ageing slides is not recommended as it may reduce signal fluorescence.

Chromoprobe Multiprobe Protocol

1. Slide preparation

i. Clean template slide

Soak the template slide for 2 minutes in 100% methanol and polish dry with a clean soft tissue.

ii. Establish correct mitotic index

It is important that the intended sample has a sufficiently high mitotic index to allow detection of chromosome abnormalities. To check the density of the sample, using a micropipette (e.g. a Gilson P10 or P20) pipette 4µl of the cell suspension onto one of the areas of the spare template slide and allow to air dry. The small volume used means that you usually have to gently touch the slide with the pipette tip to transfer the suspension. Examine by phase contrast microscopy.

If the cell density is too high, dilute the suspension with fresh fixative.

If the mitotic index is too low, spin down the fixed cell suspension at 160xg for 10 minutes. Note the volume of supernatant, remove, and re-suspend the cell pellet in a smaller volume of fresh fixative.

If cell sample concentration has been altered, spot 4 μ l of the concentrated sample onto another square of your test slide and re-examine by phase contrast microscopy.


Please Note: 50 μ l is the minimum volume required for the OctoChrome device protocol.

iii. Quality control of samples

Samples should be examined for cytoplasm since this will interfere with the *in situ* protocol. If the chromosomes appear to be enclosed by a granular material when examined under phase contrast microscopy, then this will compromise results. One method for reducing cytoplasm is to spot 4 μ l of your sample onto the template slide and watch the fixative as it spreads out. In the normal situation, the fixative will spread to maximum, recede and then evaporate. To clean up any cytoplasm we have found that effective results are achieved if a fresh drop of 3:1 fixative is allowed to fall onto the spot at the point when the spreading fixative has reached its maximum. Allow the drop of fixative to evaporate and re-examine the spot.

iv. Spotting of slide

Pipette 4 μ l of cell suspension onto all 8 areas of the template slide in a sequence of alternating squares as shown below. This will prevent the cell spreads from interfering with each other.

Once the first group of drops has air-dried, spot the remaining squares with 4 μ l drops in the same manner. After the slide has dried, examination of the slide under phase contrast will reveal whether any squares have been missed.

If spots have been missed, or squares have too few cells, simply spot those squares again : it is not necessary to re-spot a new slide.

If upon examination of slide, a square has insufficient cells/metaphases, spotting further drop(s) of suspension can be added to increase the cell density.

Please note: If the metaphase cells appear overspread then clean the template slide thoroughly in methanol and re-spot allowing every spot to dry before proceeding to the next.

2. Preparation of OctoChrome device and template slide

- i. Ensure that the Chromoprobe Multiprobe Hybridisation Chamber is in the 37°C water bath and allow to equilibrate to 37°C (+/- 1°C). (This may take up to an hour if the water bath has been switched on from cold).
- ii. Mix the hybridisation solution by repeated pipetting and pre-warm a 25µl aliquot per OctoChrome device to 37°C. Also pre-warm each OctoChrome device to 37°C by placing the device **label side down**.
Do not touch the raised boss surfaces of the OctoChrome device.
- iii. Wash template slides containing fixed samples in 2 x SSC for 2 minutes at room temperature (20 - 25°C).
- iv. Whilst the OctoChrome device is still at 37°C, dehydrate template slides containing fixed samples through an ethanol series (2 minutes each in 70%, 85% and absolute ethanol), dry and place at 37°C to warm.
- v. Add 2µl of pre-warmed hybridisation solution to each of the eight areas on the pre-warmed OctoChrome device using a P10 micropipette while it remains at 37°C.

3. Positioning of template slide over the OctoChrome device.

- i. Carefully invert the template slide over the OctoChrome device such that the number 1, which is now upside down, is located over the top right hand area of the OctoChrome device (Figure 1).

To help locate square 1 (Chromosomes 3, 15 & 17), its position on the device has been marked in Orange.


Figure 1. Location of whole chromosome painting probes on OctoChrome device squares.

- ii. Make sure that the template slide is carefully aligned with the matching areas on the OctoChrome device. Carefully lower the slide over the OctoChrome device so that the drops of hybridisation solution make contact with the slide. Apply gentle, even pressure to ensure that the hybridisation solution is spread to the edges of each of the raised areas on the OctoChrome device.
- iii. Lift the slide/OctoChrome **carefully** holding the frosted end of the glass slide and invert so that the slide is underneath the OctoChrome device. Make sure the device does not smear across the template slide as this could cause cross-contamination of the probes.

4. Instructions for use of the Cytocell slide surface thermometer.

The temperature of the hotplate should be checked with the Cytocell slide surface thermometer before proceeding to denaturation.

This thermometer is a liquid crystal device and although reversible, it must be treated with care to ensure a reasonable life span. The thermometer must only be used to check the temperature of a hotplate, it must not be used to monitor the hotplate performance over time.

To use the thermometer properly, place it onto the surface of the hotplate and wait until the different segments stop changing colour. The correct temperature is indicated by a pale green / gold colour. When the segments appear granular and the colours no longer appear uniform and regular, the thermometer should be discarded as it is exhausted. The life span of each thermometer should, however, easily be sufficient for a ten-device OctoChrome kit.

5. Denaturation

A PCR thermal cycler heating block is NOT suitable for use in place of solid bed hotplate for this procedure.

Transfer the slide/OctoChrome device to the hotplate taking particular care to hold it level. (Ensure the sample slide is in good contact with the hotplate). Denature on the hotplate at 75°C (+/- 1°C) **for 5 minutes**.

6. Hybridisation

Place slide/OctoChrome device in the Chromoprobe Multiprobe Hybridisation Chamber supplied, replace the lid and float the chamber in the 37°C (+/- 1°C) water bath (non - stirring) overnight.

Please note: **Do not seal the lid on the hybridisation chamber.**
 Do not place a lid on the water bath.
 Do not hybridise in an incubator.

Please ensure that the hybridisation chamber is completely dry. (ie. no water or damp tissue inside the chamber).

The humidity inside the chamber is vital for optimal hybridisation. The correct levels will be achieved following those steps.

7. Post-hybridisation stringent washes

i. Preparation of stringent wash solutions

1. Solution 1: Prepare a coplin/Hellendahl jar containing 0.4 x SSC. Place in a water bath and allow to reach 72°C (+/- 1°C) adjust pH to 7.0.
2. Solution 2: Prepare a coplin/Hellendahl jar containing 2 x SSC and 0.05% Tween 20. Allow standing at room temperature.

Check the temperature and pH of the solutions in coplin jars and adjust if necessary. The pH should be 7.0 when at the correct temperature.

ii. Stringent wash steps

1. Remove the OctoChrome device carefully from the slide and place the slide in Solution 1 for 2 minutes. (The OctoChrome device cannot be re-used)
2. Place the slide into Solution 2 for 30 seconds.

Avoid processing more than two OctoChrome slides through the stringency washes at any one time.

8. Mounting and visualisation of results

The slide can be either counterstained in PI or DAPI. PI will allow simultaneous visualisation of the three fluorophores using the recommended triple filter. Distinguishing the PI from the red fluorophore will not be optimal with the triple filter or the single Texas Red filter. Image analysis capturing of each fluorophore with the specific single filters (see Optimal Microscope set up) is possible, although the PI may interfere with the Texas Red channel.

With DAPI as the counterstain the blue fluorophore may appear lost when using the triple filter (because the DAPI may mask the blue fluorophore), although green and red signals will still be visible. Visual analysis and image capture are possible with the appropriate single filters. The advantage of using DAPI over PI is that the banding is superior and will allow inversion to G banding using an image analysis system.

Alternatively the slide can be stained firstly in PI and later de-stained and re-stained in DAPI. However, please note it is more difficult to remove DAPI counterstain so the initial stain should be PI.

i. Counterstaining with PI or DAPI

1. Apply 20µl of DAPI or PI to each end of the slide and apply a coverslip (24 x 50 mm)
2. Blot the slide with filter paper or tissue.

3. Leave in the dark for 10 minutes before viewing by fluorescence microscopy.
- ii. Removal of PI and re-staining with DAPI
 1. Wash slide for 5 minutes in 1 x ST.
 2. Repeat to give a total of four washes.
 3. Apply 20µl of DAPI to both ends of the slide and apply a large coverslip (24 x 50 mm)
 4. Blot the slide with filter paper.
 5. Leave in the dark for 10 minutes before viewing by fluorescence microscopy.
- iii. Certain types of microscope have slide holders which make it difficult to view the extreme ends of the slide. If this occurs then simply turn the slide through 180° which will help with the viewing of the slide.

The probes used on the Multiprobe device are directly labelled with fluorophores which are light sensitive. Results are improved when the probes are exposed to minimal amounts of light during these procedures; however, it is not necessary to work in the dark.

Stability of finished slides

FISHed slides remain analysable for up to 1 month if stored in the dark at or below room temperature.

Procedural recommendations

1. Baking or otherwise ageing of slides is not recommended as it may reduce signal fluorescence.
2. Hybridisation conditions may be adversely affected by the use of reagents other than those provided or recommended by Cytocell Technologies Ltd.
3. Failure to follow all procedures for slide preparation, denaturation and hybridisation may give unacceptable or erroneous results.
4. The use of a calibrated thermometer is strongly recommended for measuring temperatures of solutions, waterbaths, and incubators, as these temperatures are critical for optimum product performance.
5. The wash concentrations (stringency), pH and temperature are important, as low stringency can result in non-specific binding of the probe and too high stringency can result in a lack of signal.
6. Incomplete denaturation can result in lack of signal and over denaturation can also result in non-specific binding.

Expected Results

1. The acrocentric chromosome painting probes contain short arm material. This is shared between the D & G group chromosomes, so cross-hybridisation in these regions may be observed.

2. The chromosome 16 whole chromosome painting probe may show faint cross-hybridisation to the heterochromatic regions of chromosome Y.
3. The Y chromosome painting probe contains the pseudoautosomal regions common with the X chromosome. Consequently, cross-hybridisation may be observed in these regions.
4. Chromosomes 1, 5 and 19 have centromeric DNA sequences in common. Consequently cross-hybridisations in these regions may be observed between these chromosomes.
5. The interpretation of FISH results should be made in conjunction with proper controls, applying cytogenetic analytical skills and within the context of the patient's medical history and other clinical findings.
6. We recommend that the FISH probes must be used in accordance with the product label.

Customer Support

Please contact the Cytocell Sales and Marketing Department.

Probe Specifications

OctoChrome Square	Chromosome	Fluorophore
1	3	red
	15	blue
	17	green
2	8	red
	12	blue
	21*	green
3	1*	red
	16	blue
	19*	green
4	2	red
	13*	blue
	20	green
5	9	red
	11	blue
	22*	green
6	4	red
	14*	blue
	18	green
7	5*	red
	10	blue
	7	green
8	X	red
	6	blue
	Y*	green

* See the “Expected Results” paragraph for details.

References

- Collins C., Lin Kuo W., Segreaves R., Fuscoe J. Pinkel P. and Gray J.W. (1991). Construction and characterisation of plasmid libraries enriched in sequences from single human chromosomes. *Genomics* **11**: 997-1006.
- Evans J.W., Chang J.A., Giaccia A.J., Pinkel D. & Brown J.M. (1991). The use of fluorescence *in situ* hybridization combined with premature chromosome condensation for the identification of chromosome damage. *Brit. J. Cancer* **63**: 517 - 521.

- Leedham P., Lawrie M., Randhawa H., Gould C., Viegas-Pequignot E., Armstrong S. & Hulten M. (1996). Improved prenatal diagnosis by FISH Multiprobe and immunocytochemistry using a monoclonal antibody against 5-methylcytosine. *Eur. J. Hum. Genet.* **4**: Suppl 1.35
- Maher E.J., Larin Z., Smith S., Cardy D.L.N. & Southern E.M. (1995). Fluorescence *in situ* hybridisation (FISH) of twenty four different chromosome specific paints on a single microscope slide. In The Association of Clinical Cytogeneticists Abstract booklet.
- Pinkel D., Landegent J., Collins C., Fuscoe J., Segraves R., Lucas J. & Gray J. (1988). Fluorescence *in situ* hybridization with human chromosome-specific libraries: Detection of trisomy 21 and translocations of chromosome 4. *Proc. Natl. Acad. Sci. USA* **85**: 9138-9142.
- Price C. M. (1993). Fluorescence *in situ* hybridisation (review). *Blood Reviews* **7**: 127-134
- Trask B. J. (1991) Fluorescence *in situ* hybridization: applications in cytogenetics and gene mapping. *Trends in Genetics* **7**: 149-154.

Symbol Legend

REF	Catalogue Number		Consult Instructions For Use (IFU)
CONT	Contents	LOT	Lot Number
	Manufacturer		Use By (Expiry Date)
	Contains sufficient for <n> tests		Temperature Limitation (Storage Temperature; upper and lower limits)
IVD	<i>In Vitro</i> Diagnostic Medical Device		

Patents and Trademarks

Chromoprobe, Cytocell and Chromoprobe Multiprobe are registered trademarks of Cytocell Technologies Ltd.

The Chromoprobe principle is covered by international patents WO9314223, EP0623177. The design of the Multiprobe is a registered design and is also covered by a Design Patent No. 420,745.

Any cyanine dyes used in this Product are manufactured on behalf of Amersham Pharmacia Biotech Inc. under an exclusive license from Carnegie Mellon University and are covered by US Patent Number 5 268 486 and other patents pending. The Compound in this Product is manufactured by NEN Life Science Products, Inc. under US Patent Numbers 5 047 519 and 5 151 507. Use of the Product for commercial purposes is strictly forbidden without written permission from Amersham Pharmacia Biotech Inc. and NEN Life Science Products, Inc.

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Bibliographie

- Collins C., Lin Kuo W., Segraves R., Fuscoe J. Pinkel P. and Gray J.W. (1991). Construction and characterisation of plasmid libraries enriched in sequences from single human chromosomes. *Genomics* **11**: 997-1006.
- Evans J.W., Chang J.A., Giaccia A.J., Pinkel D. & Brown J.M. (1991). The use of fluorescence *in situ* hybridization combined with premature chromosome condensation for the identification of chromosome damage. *Brit. J. Cancer* **63**: 517 - 521.

- Leedham P., Lawrie M., Randhawa H., Gould C., Viegas-Pequignot E., Armstrong S. & Hulten M. (1996). Improved prenatal diagnosis by FISH Multiprobe and immunocytochemistry using a monoclonal antibody against 5-methylcytosine. *Eur. J. Hum. Genet.* 4: Suppl 1.35
- Maher E.J., Larin Z., Smith S., Cardy D.L.N. & Southern E.M. (1995). Fluorescence *in situ* hybridisation (FISH) of twenty four different chromosome specific paints on a single microscope slide. In The Association of Clinical Cytogeneticists Abstract booklet.
- Pinkel D., Landegent J., Collins C., Fuscoe J., Segraves R., Lucas J. & Gray J. (1988). Fluorescence *in situ* hybridization with human chromosome-specific libraries: Detection of trisomy 21 and translocations of chromosome 4. *Proc. Natl. Acad. Sci. USA* 85: 9138-9142.
- Price C. M. (1993). Fluorescence *in situ* hybridisation (review). *Blood Reviews* 7: 127-134
- Trask B. J. (1991) Fluorescence *in situ* hybridization: applications in cytogenetics and gene mapping. *Trends in Genetics* 7: 149-154.

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Bibliografia

- Collins C., Lin Kuo W., Segaves R., Fuscoe J. Pinkel P. and Gray J.W. (1991). Construction and characterisation of plasmid libraries enriched in sequences from single human chromosomes. *Genomics* **11**: 997-1006.
- Evans J.W., Chang J.A., Giaccia A.J., Pinkel D. & Brown J.M. (1991). The use of fluorescence *in situ* hybridization combined with premature chromosome condensation for the identification of chromosome damage. *Brit. J. Cancer* **63**: 517 - 521.

- Leedham P., Lawrie M., Randhawa H., Gould C., Viegas-Pequignot E., Armstrong S. & Hulten M. (1996). Improved prenatal diagnosis by FISH Multiprobe and immunocytochemistry using a monoclonal antibody against 5-methylcytosine. *Eur. J. Hum. Genet.* **4**: Suppl 1.35
- Maher E.J., Larin Z., Smith S., Cardy D.L.N. & Southern E.M. (1995). Fluorescence *in situ* hybridisation (FISH) of twenty four different chromosome specific paints on a single microscope slide. In The Association of Clinical Cytogeneticists Abstract booklet.
- Pinkel D., Landegent J., Collins C., Fuscoe J., Seagraves R., Lucas J. & Gray J. (1988). Fluorescence *in situ* hybridization with human chromosome-specific libraries: Detection of trisomy 21 and translocations of chromosome 4. *Proc. Natl. Acad. Sci. USA* **85**: 9138-9142.
- Price C. M. (1993). Fluorescence *in situ* hybridisation (review). *Blood Reviews* **7**: 127-134
- Trask B. J. (1991) Fluorescence *in situ* hybridization: applications in cytogenetics and gene mapping. *Trends in Genetics* **7**: 149-154.

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Literatur

- Collins C., Lin Kuo W., Segraves R., Fuscoe J. Pinkel P. and Gray J.W. (1991). Construction and characterisation of plasmid libraries enriched in sequences from single human chromosomes. *Genomics* **11**: 997-1006.
- Evans J.W., Chang J.A., Giaccia A.J., Pinkel D. & Brown J.M. (1991). The use of fluorescence *in situ* hybridization combined with premature chromosome condensation for the identification of chromosome damage. *Brit. J. Cancer* **63**: 517 - 521.
- Leedham P., Lawrie M., Randhawa H., Gould C., Viegas-Pequignot E., Armstrong S. & Hulten M. (1996). Improved prenatal diagnosis by FISH Multiprobe and

immunocytochemistry using a monoclonal antibody against 5-methylcytosine. *Eur. J. Hum. Genet.* **4**: Suppl 1.35

- Maher E.J., Larin Z., Smith S., Cardy D.L.N. & Southern E.M. (1995). Fluorescence *in situ* hybridisation (FISH) of twenty four different chromosome specific paints on a single microscope slide. In The Association of Clinical Cytogeneticists Abstract booklet.
- Pinkel D., Landegent J., Collins C., Fuscoe J., Segraves R., Lucas J. & Gray J. (1988). Fluorescence *in situ* hybridization with human chromosome-specific libraries: Detection of trisomy 21 and translocations of chromosome 4. *Proc. Natl. Acad. Sci. USA* **85**: 9138-9142.
- Price C. M. (1993). Fluorescence *in situ* hybridisation (review). *Blood Reviews* **7**: 127-134
- Trask B. J. (1991) Fluorescence *in situ* hybridization: applications in cytogenetics and gene mapping. *Trends in Genetics* **7**: 149-154.

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- Collins C., Lin Kuo W., Segreaves R., Fuscoe J., Pinkel P., and Gray J.W. (1991). Construction and characterisation of plasmid libraries enriched in sequences from single human chromosomes. *Genomics* **11**: 997-1006.
- Evans J.W., Chang J.A., Giaccia A.J., Pinkel D. & Brown J.M. (1991). The use of fluorescence *in situ* hybridization combined with premature chromosome condensation for the identification of chromosome damage. *Brit. J. Cancer* **63**: 517 - 521.
- Leedham P., Lawrie M., Randhawa H., Gould C., Viegas-Pequignot E., Armstrong S. & Hulten M. (1996). Improved prenatal diagnosis by FISH Multiprobe and immunocytochemistry using a monoclonal antibody against 5-methylcytosine. *Eur. J. Hum. Genet.* **4**: Suppl 1.35.

- Maher E.J., Larin Z., Smith S., Cardy D.L.N. & Southern E.M. (1995). Fluorescence *in situ* hybridisation (FISH) of twenty four different chromosome specific paints on a single microscope slide. In The Association of Clinical Cytogeneticists Abstract booklet.
- Pinkel D., Landegent J., Collins C., Fuscoe J., Segraves R., Lucas J. & Gray J. (1988). Fluorescence *in situ* hybridization with human chromosome-specific libraries: Detection of trisomy 21 and translocations of chromosome 4. Proc. Natl. Acad. Sci. USA **85**: 9138-9142.
- Price C. M. (1993). Fluorescence *in situ* hybridisation (review). Blood Reviews **7**: 127-134
- Trask B. J. (1991) Fluorescence *in situ* hybridization: applications in cytogenetics and gene mapping. Trends in Genetics **7**: 149-154.

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

Referencer

- Collins C., Lin Kuo W., Segraves R., Fuscoe J. Pinkel P. and Gray J.W. (1991). Construction and characterisation of plasmid libraries enriched in sequences from single human chromosomes. *Genomics* **11**: 997-1006.
- Evans J.W., Chang J.A., Giaccia A.J., Pinkel D. & Brown J.M. (1991). The use of fluorescence *in situ* hybridization combined with premature chromosome condensation for the identification of chromosome damage. *Brit. J. Cancer* **63**: 517 - 521.
- Leedham P., Lawrie M., Randhawa H., Gould C., Viegas-Pequignot E., Armstrong S. & Hulten M. (1996). Improved prenatal diagnosis by FISH Multiprobe and immunocytochemistry using a monoclonal antibody against 5-methylcytosine. *Eur. J. Hum. Genet.* **4**: Suppl 1.35
- Maher E.J., Larin Z., Smith S., Cardy D.L.N. & Southern E.M. (1995). Fluorescence *in situ* hybridisation (FISH) of twenty four different chromosome specific paints on a single microscope slide. In *The Association of Clinical Cytogeneticists Abstract booklet*.
- Pinkel D., Landegent J., Collins C., Fuscoe J., Segraves R., Lucas J. & Gray J. (1988). Fluorescence *in situ* hybridization with human chromosome-specific libraries: Detection of trisomy 21 and translocations of chromosome 4. *Proc. Natl. Acad. Sci. USA* **85**: 9138-9142.
- Price C. M. (1993). Fluorescence *in situ* hybridisation (review). *Blood Reviews* **7**: 127-134
- Trask B. J. (1991) Fluorescence *in situ* hybridization: applications in cytogenetics and gene mapping. *Trends in Genetics* **7**: 149-154.

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Chromoprobe, Cytocell og Chromoprobe Multiprobe er registrerede varemærker for Cytocell Technologies Ltd.