


## SuperLight™ Luciferase Reporter Gene Assay Kits

### Bioluminescent Assay for Promoter Regulated Luciferase Expression

#### DESCRIPTION

The SuperLight™ Luciferase Reporter Gene Assay is based on the quantitation of luciferase expression in mammalian, yeast or E. coli cells, using luciferin and ATP as substrates. The reaction results in light production which can be conveniently measured on a luminometer.

##### Reaction scheme:



This bioluminescent reporter gene assay is extremely sensitive and is especially suitable for quantifying luciferase expression in recombinant cells. This ultra-sensitive, homogeneous cell-based assay only requires adding a single reagent to the cells and measuring the light intensity after a short incubation step (2 minutes). Assays can be performed in tubes, cuvettes or multi-well plates. All kit components are compatible with culture media and with all liquid handling systems. With an extended luminescence emission kinetics (half-life 40 min), the SuperLight™ luciferase assays are especially suitable for high-throughput screening in 96-well, 384-well and 1536-well plates. In addition, the reagent provided in the kits has been formulated for maximum sensitivity, reproducibility and long shelf-life. Applications for this kit include gene regulation studies and high-throughput screening of gene modulators.

#### KEY FEATURES

**High sensitivity and wide detection range:** detection of as little of 2 fg luciferase and as few as 4 cells. Plus, the emitted light is linear over seven orders of magnitude.

**Compatible with routine laboratory and HTS formats:** assays can be performed in tubes or microplates, on LJL Analyst, Berthold Luminometer, Top-Count, MicroBeta counters, chemiluminescent image plate readers (CLIPR/LeadSeeker). Assay reagents compatible with all liquid handling systems.

**Fast and convenient:** homogeneous “mix-and-measure” assay allows detection of luciferase levels within 10 minutes. The optimally combined reagent system allows a single addition step, and simultaneous cell lysis and detection.

**Robust and amenable to HTS:** Z' factors of 0.6 to 0.8 are observed in 96-well and 384-well plates. Can be readily automated on HTS liquid handling systems.

#### APPLICATIONS

**Gene Regulation:** gene expression level, characterization of promoter activity, modulation of gene expression by receptors, transcription factors and small molecules.

**Drug Discovery:** high-throughput screen for gene modulators.

#### KIT CONTENTS

Catalog #	Size (assays)	Reagent	Assay Buffer
SLLU-200	200	solid	20 mL
SLLU-500	500	solid	50 mL
SLLU-01K	1,000	solid	100 mL
SLLU-HTS	> 5,000	solid	customized

**Storage conditions.** Store the Reagent in the provided amber tube at -20°C and the Assay Buffer at 2-8°C. Shelf life: 12 month.

**This protocol** can be downloaded online at [www.bioassaysys.com](http://www.bioassaysys.com).

**Precautions:** reagents are for research use only. Normal precautions for laboratory reagents should be exercised while using the reagents. Please refer to Material Safety Data Sheet for detailed information.

#### PROCEDURES

The SuperLight™ Luciferase Reporter Gene Assay is based on the bioluminescence generated during the luciferin/luciferase reaction. The reconstituted reagent has been optimized to combine cell lysis and detection into one single step. Phenol red in culture media does not interfere in this assay. All data in the Technical Notes were obtained in media containing phenol red.

**Important:** fresh reconstitution of the Reagent in Assay Buffer is recommended, although the reconstituted Reagent may be stable for up to 4 weeks when stored at -20 °C.

#### Procedure using 96-well plate:

1. Plate and culture cells (80 µL) in white opaque 96-well tissue culture plates. Typical culture medium contains DMEM, 10% fetal bovine serum and antibiotics (penicillin/ streptomycin, gentamycin, etc). Amino acids and other nutrients can be added to the culture medium. Assays can be performed on either adherent cells or cells in suspension. The cells can be either stably or transiently transfected with the luciferase gene. Culture volume can vary from 50 to 100 µL, although 80 µL is used in this protocol. Blank control wells containing no cells should also be prepared.

2. Add test compounds and controls to cells. Mix well and incubate for the cells desired period of time. Incubation time for gene regulation studies can be from several hours up to 3 days. It is recommended that assays be run in duplicate or triplicate. A volume of 20 µL compounds in PBS or culture medium is recommended.

3. Reconstitute the Reagent. First equilibrate the Reagent and Assay Buffer to room temperature. Then simply combine the Assay Buffer and Reagent by pipetting a small volume (e.g. 1 mL) buffer to the Reagent tube. Vortex briefly and pipet the reconstituted solution to the Assay Buffer bottle. Repeat this step to transfer all Reagent to the Assay Buffer bottle. Mix by inversion until the Reagent is thoroughly dissolved. After this is done, mark the bottle label as Reconstituted Reagent.

4. Add 100 µL (per 80 µL of cell culture) of the reconstituted Reagent to each well and mix well with the cells. Incubate for 2 minutes at room temperature. The volume of the reagent can be adjusted depending on the volume of cell culture.

5. Measure luminescence on a luminometer. The integration time can be 1 sec to 2 min depending on the luciferase expression level and instrument sensitivity. For most luminometers (Berthold Luminometer, LJL Analyst, Top Count, MicroBeta Counters, CLIPR and LeadSeeker), integration 1 to 5 sec is appropriate.

#### Procedure using 384-well plate:

1. Plate and culture cells (25 µL) in white opaque 384-well tissue culture plates. Culture volume can vary from 20 to 50 µL, although 25 µL is used in this protocol. Set up blank control wells containing no cells.

2. Add test compounds and controls to cells. Mix well and incubate for the cells desired period of time. A volume of 5 µL compounds in PBS or culture medium is recommended.

3. Reconstitute the Reagent using the same procedure as for the 96-well assay.

4. Add 30 µL (per 25 µL of cell culture) of the reconstituted Reagent per well and mix well with the cells. Incubate for 2 minutes at room temperature. The volume of the reagent can be adjusted depending on the volume of cell culture.

5. Measure luminescence on a luminometer. The integration time can be 1 sec to 2 min depending on the luciferase expression level and instrument sensitivity. For most luminometers (Berthold Luminometer, LJL Analyst, Top Count, MicroBeta Counters, CLIPR and LeadSeeker), integration 1 to 5 sec is appropriate.

#### GENERAL CONSIDERATIONS

**Incubation time.** Both the luciferin/luciferase reaction and cell lysis are fast, so incubation for 2 to 10 minutes following reagent addition is generally enough for mammalian cells (e.g. HEK293, CHO).

**Cell number.** The optimized reporter gene assay reagent is very sensitive to luciferase (detection limit 2 fg) and exhibits linearity over seven orders of magnitude. As few as 4 cells can be determined and a linear response is still observed with as many as 80,000 cells per 96-well. For assay optimization, it is recommended that the optimal number

## SuperLight™ Luciferase Reporter Gene Assay Kits

### Bioluminescent Assay for Promoter Regulated Luciferase Expression

of cells per well be determined by serial dilution of cells. Cells can be adherent or in suspension cultures.

**Cell lysis and mixing.** For the sake of convenience, the addition of 1 volume of reconstituted reagent to 1 volume of cells allows a sufficient mixing. No additional mixing is required since the specially formulated buffer instantly lyses mammalian cells.

#### DATA ANALYSIS

The light intensity (RLU) is directly proportional to the luciferase concentration. For dose-response studies, the data are plotted against compound concentration and the EC<sub>50</sub> for gene up-regulator compound and IC<sub>50</sub> for a gene down regulator compound can be determined by non-linear regression analysis using Prism or other data analysis tools.

#### LITERATURE

1. Zhao L, Haslam DB (2005). A quantitative and highly sensitive luciferase-based assay for bacterial toxins that inhibit protein synthesis. *J Med Microbiol.* 54:1023-30. *A most recent publication of use of our reagent.*

#### Gene Regulation (receptors, small molecules, transcription factors)

- Imaki H et al (2003). Cell cycle-dependent regulation of the Skp2 promoter by GA-binding protein. *Cancer Res.* 63(15):4607-13.
- Nabors LB et al (2003). Tumor necrosis factor alpha induces angiogenic factor up-regulation in malignant glioma cells: a role for RNA stabilization and HuR. *Cancer Res.* 63(14):4181-7.
- Keiss HP et al (2003). Garlic (*Allium sativum* L.) modulates cytokine expression in lipopolysaccharide-activated human blood thereby inhibiting NF-kappaB activity. *J Nutr.* 133(7):2171-5.

#### Characterization of promoter activity

- Sharina IG et al (2003). CCAAT-binding factor regulates expression of the beta1 subunit of soluble guanylyl cyclase gene in the BE2 human neuroblastoma cell line. *Proc Natl Acad Sci U S A.* 100(20):11523-11528.
- Cramer SD et al (2003). Association between genetic polymorphisms in the prostate-specific antigen gene promoter and serum prostate-specific antigen levels. *J Natl Cancer Inst.* 95(14):1044-53.

#### Modulation of gene expression by small molecules

- Kim SH et al (2003). Luteolin inhibits the nuclear factor-kappa B transcriptional activity in Rat-1 fibroblasts. *Biochem Pharmacol.* 66:955-63.

#### High-throughput screening for gene modulators Reviews:

- Phippard D, Manning AM (2003). Screening for inhibitors of transcription factors using luciferase reporter gene expression in transfected cells. *Methods Mol Biol.* 225:19-23.
- Naylor LH (1999). Reporter gene technology: the future looks bright. *Biochem Pharmacol.* 58:749-57.

#### cAMP response element (CRE), NFAT and SRE dependent pathways:

- Chen G et al (1999). Constitutive receptor systems for drug discovery. *J Pharmacol Toxicol Methods.* 42:199-206.
- Baldari CT et al (1998). NF-AT-luciferase reporter T cell lines as tools to screen immunosuppressive drugs. *Biologicals.* 26:1-5.
- Choi SE, Choi EY, Kim PH, Kim JH. Involvement of protein kinase C and rho GTPase in the nuclear signalling pathway by transforming growth factor-beta1 in rat-2 fibroblast cells. *Cell Signal.* 1999;11(1):71-6.
- Rapisarda A et al (2002). Identification of small molecule inhibitors of hypoxia-inducible factor 1 transcriptional activation pathway. *Cancer Res.* 62:4316-24.

#### TECHNICAL NOTES

The SuperLight™ Luciferase Reporter Gene Assay Kit has been specially optimized and formulated to provide a sensitive, convenient and robust assay for gene expression and regulation studies in mammalian cells. Key

features of the kits are as follows:

**High sensitivity and wide detection range:** detection of as little of 2 fg luciferase and as few as 4 cells. Plus, the emitted light is linear over seven orders of magnitude.

**Compatible with routine laboratory and HTS formats:** assays can be performed in tubes or microplates, on LJL Analyst, Berthold Luminometer, Top-Count, MicroBeta counters, chemiluminescent image plate readers (CLIPR/LeadSeeker). Assay reagents compatible with all liquid handling systems.

**Fast and convenient:** homogeneous "mix-and-measure" assay allows detection of luciferase levels within 10 minutes. The optimally combined reagent system allows a single addition step and simultaneous cell lysis and detection.

**Robust and amenable to HTS.** Z' factors of 0.7 to 0.9 are routinely observed in 96-well and 384-well plates. Can be readily automated on HTS liquid handling systems.




Figure 1. Linearity of the SuperLight™ Luciferase Assay in 384-well plate. The detection limit estimated from the blank controls was 2 fg. Light intensity was linear from 2 fg up to 46 ng.




Figure 2. Linear relationship between emitted light and number of HEK293 cells transiently transfected with a CRE-luciferase reporter construct in a 384-well plate assay. Detection limit: 4 cells.




Figure 3. Up-regulation of CRE-dependent luciferase expression by prostaglandin E2 (PGE<sub>2</sub>) and adenylyl cyclase activator forskolin in HEK293 cells transiently transfected with a CRE-luciferase construct. EC<sub>50</sub> = 0.15 μM for PGE<sub>2</sub> and 3.5 μM for forskolin.