

# InVitroCare® IVF MEDIA SERIES

## HTF HEPES

For **egg retrieval, follicle flushing**, ICSI procedures and cryopreservation

Formulated with EDTA and a stabilized form of glutamine  
18 month shelf life at 2-8°C  
Cat # 2002 - 100ml or Cat # 2002-5 - 500ml

## HTF

For use as a **fertilization media** and for early embryo culture through 2PN stage

Formulated with EDTA and a stabilized form of glutamine  
120 day shelf life at 2-8°C  
Cat # 2001 - 100ml

## IVC-ONE™

For **early embryo development** from Day 0 through compaction

**Glucose and Phosphate free**  
Formulated with EDTA, Citrate, Glutathione, and a stabilized form of glutamine  
120 day shelf life at 2-8°C  
Cat # 2006 - 100ml

## IVC-TWO™

For use in **Day 3 embryo transfers** and ICSI procedures

**Low Glucose and Phosphate free**  
Formulated with EDTA, Citrate, Glutathione, and a stabilized form of glutamine  
120 day shelf life at 2-8°C  
Cat # 2008 - 100ml

## IVC-THREE™

For **extended embryo culture** and transfer through expanded **Blastocyst (Day 5)**

**Phosphate free and elevated levels of Glucose**  
Supplemented with amino acids and vitamins  
Formulated with Glutathione  
120 day shelf life at 2-8°C  
Cat # 2007 - 100ml