

## EUROCLONE MOLECULAR DIAGNOSTICS INFECTIOUS DISEASES PRODUCT LINE

### Infectious Diseases Respiratory panel

Code	Description	Format	Shelf Life	Availability	Technology	Certification	Flexible Run Size	Target
EBR001032	Duplic $\alpha$ RealTime Mycobacterium tuberculosis	32 tests	12 months	available	Real time PCR	CE-IVD	4 runs_ 6 samples and 2 CTRLs each one	16-23s rRNA gene
EBR010032	Duplic $\alpha$ RealTime Mycobacterium Avium complex	32 tests	12 months	available	Real time PCR	CE-IVD	4 runs_ 6 samples and 2 CTRLs each one	16-23s rRNA gene
EBR008032	Duplic $\alpha$ RealTime Adenovirus	32 tests	12 months	available	Real time PCR	CE-IVD	4 runs_ 6 samples and 2 CTRLs each one	Exon gene
EBR007032	Duplic $\alpha$ RealTime Legionella pneumophila	32 tests	12 months	available	Real time PCR	CE-IVD	4 runs_ 6 samples and 2 CTRLs each one	23s ribosomal RNA gene
EBR011032	Duplic $\alpha$ RealTime Legionella ssp	32 tests	12 months	available	Real time PCR	CE-IVD	4 runs_ 6 samples and 2 CTRLs each one	16s ribosomal RNA gene
EBR006032	Duplic $\alpha$ RealTime Chlamydia pneumoniae	32 tests	12 months	available	Real time PCR	CE-IVD pending	4 runs_ 6 samples and 2 CTRLs each one	MOMP gene
EBR012032	Duplic $\alpha$ RealTime Mycoplasma pneumoniae	32 tests	12 months	available	Real time PCR	CE-IVD	4 runs_ 6 samples and 2 CTRLs each one	Mp1842 cytheadhesion P1 gene

### Infectious Diseases Sex. Transm. panel

Code	Description	Format	Shelf Life	Availability	Technology	Certification	Flexible Run Size	Target
EBR005032	Duplic $\alpha$ RealTime Chlamydia trachomatis	32 tests	12 months	available	Real time PCR	CE-IVD pending	4 runs_ 6 samples and 2 CTRLs each one	16S rRNA gene
EBR009032	Duplic $\alpha$ RealTime Thricomonas vaginalis	32 tests	12 months	available	Real time PCR	CE-IVD	4 runs_ 6 samples and 2 CTRLs each one	Beta-tubulin (btub3) gene
EBR013032	Duplic $\alpha$ RealTime Neisseria gonorrhoeae	32 tests	12 months	available	Real time PCR	CE-IVD	4 runs_ 6 samples and 2 CTRLs each one	Cryptic plasmid pJD1
EBR014032	Duplic $\alpha$ RealTime Mycoplasma genitalium	32 tests	12 months	available	Real time PCR	CE-IVD	4 runs_ 6 samples and 2 CTRLs each one	Repetitive sequence element mgp-r2
EBR015032	Duplic $\alpha$ RealTime Mycoplasma hominis	32 tests	12 months	available	Real time PCR	CE-IVD	4 runs_ 6 samples and 2 CTRLs each one	16s ribosomal RNA gene
EBR016032	Duplic $\alpha$ RealTime Ureaplasma urealyticum	32 tests	12 months	available	Real time PCR	CE-IVD	4 runs_ 6 samples and 2 CTRLs each one	Urease complex component



EUROCLONE S.p.A.

Via Lombardia, 12 - 27010 Siziano (PV) Italy

Phone: +39.02.38195378 - Fax: +39.02.38195248

e-mail: info@euroclone.net - www.euroclone.net

Euroclone S.p.A. has a Quality System certified Dasa Raegister spa in compliance with UNI EN ISO 9001 (2000)

